

**PROCEEDINGS
OF THE
AUSTRALIAN
PHYSIOLOGICAL
SOCIETY**



PROCEEDINGS OF THE AUSTRALIAN PHYSIOLOGICAL SOCIETY

Volume 36

2005



Correspondence

Correspondence concerning literary and technical aspects of the journal should be addressed to the Editor. Correspondence concerning subscriptions to the journal should be addressed to the Treasurer and that involving all other Society business should be addressed to the National Secretary.

**AUSTRALIAN AUSTRALIAN PHYSIOLOGICAL
SOCIETY INC.**

President: D.J. Adams (2004-)
School of Biomedical Sciences
University of Queensland
St. Lucia
Queensland 4072

National Secretary: D.A. Saint (2002-)
Molecular and Biomedical Sciences
University of Adelaide
Adelaide SA 5005

Treasurer: C.B. Neylon (2002-)
Department of Anatomy & Cell Biology
University of Melbourne
VIC 3010

Editor: D.F. Davey (2002-)
D'Entrecasteaux
378 Manuka Road
Kettering,
Tasmania 7155

Webmaster: D.R. Laver (2005-)
School of Biomedical Sciences
University of Newcastle
NSW 2308

Elected Members of Council (Year of retirement in parentheses)

C. Chen	(2006)	Prince Henry's Institute of Medical Research
L. Delbridge	(2005)	University of Melbourne
A. Mihailidou	(2006)	Royal North Shore Hospital
D. Van Helden	(2006)	University of Newcastle
E. Porrello	(2007)	University of Melbourne (Student representative)
G. Lamb	(2008)	La Trobe University
G. Lynch	(2008)	University of Melbourne
P. Poronnik	(2008)	University of Queensland

Co-opted Members of Council

T. Lewis (2006) University of NSW (Newsletter)

© 2005 Australian Physiological Society

Sustaining Members

ADInstruments

6/4 Gladstone Rd
CASTLE HILL, NSW 2154
Phone: (02) 9899 5455
Fax: (02) 9899 5847
Email: hlalevski@adi.com.au
Web: www.adInstruments.com
Helen Lalevski

Blackwell Science Pty Ltd

PO Box 378
Carlton South, Victoria 3053
Phone: (03) 9347 0300
Fax: (03) 9347 5001
Email: suneel.jethani@blackwellpublishingasia.com
Web: www.blackwellpublishingasia.com/
Mr Suneel Jethani

Olympus Australia

PO Box 985
Mount Waverley VIC 3149
Phone 1300 132 992
Fax (03) 9543 1350
Email:
customerservice@olympus.com.au
Web: www.olympus.com.au
Violetta Mironova (Sales Manager)

Sapphire Bioscience Pty Ltd

Suite 1, 134 Redfern Street
Redfern NSW 2016
Australia
Phone: +61 2 9698 2022
Fax: +61 2 9698 1022
Email: sales@sapphirebioscience.com
Web: www.sapphirebioscience.com
Sue Goodman

SDR Clinical Technology

213 Eastern Valley Way
Middle Cove, NSW 2068
Phone: (02) 9958 2688
Fax: (02) 9958 2655
Email: sdr@sdr.com.au
Web: www.sdr.com.au
Dr Peter Kenny
Dr Roger Lainson

Taylor-Wharton (Australia) Pty Ltd

Unit 1, 882 Lesley Drive
Albury, NSW 2640
Phone: (02) 6040 2533
(Toll free - 1800 804 037)
Fax: (02) 6040 2510
Email: bsmitstw@ozemail.com.au
Mr Bill Smits (General Manager)

PROCEEDINGS OF THE
AUSTRALIAN PHYSIOLOGICAL SOCIETY

Volume 36

2005

CONTENTS

AuPS/AHMRC Sydney 2004 Meeting

British Physiology Society Lecture

Emerging functions of 10 types of transient receptor potential (TRP) cationic channel in vascular smooth muscle

D.J. Beech 1

Symposium - Environmental and genetic influences on respiratory health

Influence of air pollution on respiratory health during perinatal development

K.E. Pinkerton and J.P. Joad 9

Role of the physico-chemical environment in lung development

S.B. Hooper and M.J. Wallace 15

The respiratory consequences of preterm birth

T.J.M. Moss 23

Environmental factors and gene-environment interactions in the aetiology of asthma

G.B. Marks 29

Symposium - Cardiac hypertrophy and failure - cellular and signalling defects

Cardiac hypertrophy, substrate utilisation and metabolic remodelling: cause or effect?

R.H. Ritchie and L.M.D. Delbridge 35

Convergence of glucose- and fatty acid-induced abnormal myocardial excitation-contraction coupling and insulin signalling

A.J. Davidoff 45

Reactive oxygen species in cardiac signalling - from mitochondria to plasma membrane ion channels

L.C. Hool 55

Genes, calcium and modifying factors in hypertrophic cardiomyopathy

T. Tsoutsman, L. Lam and C. Semsarian 63

Symposium - Integrative aspects of human muscle performance

Effects of exercise-induced arterial hypoxemia on limb muscle fatigue and performance

L.M. Romer and J.A. Dempsey 71

Exercise and skeletal muscle GLUT4 expression: molecular mechanisms

S.L. McGee and M. Hargreaves 77

Evidence for a supraspinal contribution to human muscle fatigue

J.L. Taylor, G. Todd and S.C. Gandevia 83

Symposium - Stretch Activated Ion Channels

<i>Gene expression of stretch activated channels and mechano-electric feedback in the heart</i> D. Kelly, L. Mackenzie, P. Hunter, B. Smaill and D.A. Saint.....	91
<i>Mechanosensitive ion channels in skeletal muscle: A link in the membrane pathology of muscular dystrophy</i> J.B. Lansman and A. Franco-Obregón.....	101
<i>Muscle damage in mdx (dystrophic) mice: The role of calcium and reactive oxygen species</i> N.P. Whitehead, E.W. Yeung and D.G. Allen	111

AuPS/ASB Canberra 2005 Meeting - Abstracts

Meeting Opening

<i>How does the membrane electric field gate an ion channel open and close?</i> Francisco Bezanilla	1P
--	----

Symposium 1: Ion Channel Gating

<i>The domains in the Na channel have specific functions</i> Francisco Bezanilla and Baron Chanda	2P
<i>Modulation of potassium channel conformation and function by permeating ions</i> S.J. Korn.....	3P
<i>The gating of mechanosensitive ion channels</i> B. Martinac.....	4P
<i>The mechanism of fast gating in ClC chloride channels</i> B. Corry, D. Bisset and S.H. Chung.....	5P
<i>Conformational changes associated with glycine receptor activation</i> J.W. Lynch and R.L. Hawthorne	6P

Symposium 2: Functional Roles of Potassium Channels in the Vasculature

<i>Endothelium-dependent vasodilatation: Fundamental role of SK_{Ca} and IK_{Ca} potassium channels</i> C.J. Garland.....	7P
<i>Functional effects of vascular K_{IR} channels</i> C.G. Sobey	8P
<i>K_v as a target for nitroxyl anion (NO⁻)-mediated vasodilatation</i> B.K. Kemp-Harper and J.L. Favaloro	9P
<i>The smooth muscle BK_{Ca} potassium channel and its interaction with arteriolar myogenic tone</i> T.V. Murphy, Y.T. Hwang, H. Ding, N. Kotecha and M.A. Hill	10P
<i>Potassium channels in vascular dysfunction</i> C.R. Triggle, A. Ellis, L. Ceroni, W. Wiehler and H. Ding	11P

Free Communications 1: Exercise Physiology

- Sudomotor responses during isometric exercise appear to be intensity- and muscle mass-dependent*
C.J. Gordon, C.D. Haley, J.N. Caldwell and N.A.S. Taylor 12P
- Hydration indices in exertional heat stress*
A.T. Garrett, N.G. Goossens, N.J. Rehrer, M.J. Patterson and J.D. Cotter 13P
- Plasma ammonia responses during heavy-intensity constant-load cycling in young and older individuals*
S. Sabapathy, D.A. Schneider and N.R. Morris 14P
- Abnormal muscle Na^+ , K^+ -pumps, plasma K^+ , and exercise limitation in renal failure patients*
A.C. Petersen, M.J. Leikis, K.T. Murphy, J.A. Leppik, X. Gong, A.B. Kent, L.P. McMahon and M.J. McKenna 15P
- The effect of eccentric exercise on plasma K^+ regulation and skeletal muscle Na^+ , K^+ -ATPase activity and content*
J.A. Bennie, C.A. Goodman, M.J. Leikis and M.J. McKenna 16P
- N-acetylcysteine infusion enhances skeletal muscle Na^+ , K^+ -ATPase activity and plasma K^+ regulation, and delays fatigue, during prolonged submaximal exercise in well-trained individuals*
C.A. Goodman, I. Medved, M.J. Brown, A.R. Bjorksten, K.T. Murphy, A.C. Petersen, S. Sostaric, X. Gong and M.J. McKenna 17P

Free Communications 2: Ion channels

- Investigating the mechanism of proton transfer through the bacterial ClC transporter*
M. O'Mara, J. Yin, M. Hoyles and S.H. Chung 18P
- An electrostatic basis for valence selectivity in cationic channels*
T. Vora, B. Corry and S.H. Chung 19P
- A current source and a cation conductance are components of an electrical circuit connected across the plasma membrane of the malaria parasite Plasmodium falciparum*
R.J.W. Allen, K.J. Saliba and K. Kirk 20P
- Role of protein flexibility in gramicidin A channel permeability*
Turgut Bastug 21P
- Ca^{2+} influx through store-operated Ca^{2+} channel in mouse sinoatrial node*
Y.K. Ju, H. Chaulet, R.M. Graham and D.G. Allen 22P
- A hydrogen peroxide insult causes sustained alteration in the sensitivity of the L-type Ca^{2+} channel to β -adrenergic receptor stimulation in ventricular myocytes*
L.C. Hool, H.M. Viola and P.G. Arthur 23P

Plenary Lecture - Nigel Unwin

- Structure and gating mechanism of the nicotinic acetylcholine receptor*
N. Unwin 24P

Free Communications 3: Ligand-Gated Ion Channels

- Analysis of a GABA_A $\gamma 2$ (R43Q) knock-in mouse model of familial epilepsy*
S. Petrou, H. Tan, P. Davies and S. Murphy 25P

<i>The receptor-associated protein, rapsyn, and regulation of postsynaptic acetylcholine receptor packing density and turnover at the neuromuscular synapse</i> W.D. Phillips and O.L. Gervásio.....	26P
<i>Functional consequences of clustering GABA_A receptors</i> M.L. Tierney, T.L. Luu, A.E. Everitt and P.W. Gage.....	27P
<i>The charge of the P-loop glutamate controls cation-anion selectivity in CNG channels</i> W. Qu, A.J. Moorhouse, M. Chandra, K.D. Pierce, T.M. Lewis and P.H. Barry.....	28P
<i>Structure and dynamics of the periplasmic loop of the MscL mechanosensitive channel studied by electron paramagnetic resonance spectroscopy</i> G. Meyer, E. Perozo and B. Martinac.....	29P
Free Communications 4: Skeletal Muscle Regulation: From Molecular Mechanism to Physiology	
<i>Calcium-phosphate precipitation in the sarcoplasmic reticulum reduces action potential-mediated Ca²⁺ release in mammalian skeletal muscle</i> T.L. Dutka, L. Cole and G.D. Lamb.....	30P
<i>Digoxin effects on muscle strength, fatigue and K⁺ fluxes during exercise in healthy young adults</i> M.J. McKenna, S. Sostaric, M.J. Brown, C.A. Goodman, X. Gong, A.C. Petersen, J. Aw, J.A. Leppik, C.H. Steward, S.F. Fraser, R.J. Snow and H. Krum.....	31P
<i>The peak tetanic force-[K⁺]_o relationship in mouse fast- and slow-twitch muscle: modulation with [Na⁺]_o or [Ca²⁺]_o</i> S.P. Cairns.....	32P
<i>The effect of dithiothreitol (DTT) application on isolated mouse muscle fatigued at 37°C</i> T.R. Moopanar and D.G. Allen.....	33P
<i>Cytoplasmic ATP-sensing CBS domains regulate gating of skeletal muscle ClC-1 chloride channels</i> B. Bennetts, G.Y. Rychkov, H-L. Ng, C.J. Morton, D. Stapleton, M.W. Paarker and B.A. Cromer.....	34P
<i>Modelling diffusive O₂ supply to isolated muscle preparations</i> C.J. Barclay.....	35P
Free Communications presented as Posters	
Ion channel structure and function	
<i>Effects of gadolinium and static magnetic fields on MscL channel activity</i> E. Petrov, Z.-W. Liu and B. Martinac.....	36P
<i>Conformational changes involved in MscL channel gating measured using FRET spectroscopy</i> B. Corry, P. Rigby and B. Martinac.....	37P
<i>C-terminal charged cluster of the mechanosensitive channel MscL, RKKEE, functions as a pH sensor</i> A. Kloda and B. Martinac.....	38P
<i>The effects of eriochrome cyanine R on the mechanosensitive channels of E. coli</i> T. Nguyen, B. Clare, L.C. Hool and B. Martinac.....	39P
<i>Mutations within the selectivity filter of the NMDA receptor channel influence voltage-dependent block by extracellular 5-Hydroxytryptamine</i> A. Kloda and D.J. Adams.....	40P

<i>Ion selectivity of glycine receptors with mutations of charged residues in the intracellular portals</i> T.M. Lewis, Sugiharto, J.A. Peters, J.J. Lambert, P.H. Barry and A.J. Moorhouse.....	41P
<i>A molecular determinant of tropisetron inhibition of the glycine receptor Cl⁻ channel</i> Z. Yang, A.D. Ney and J.W. Lynch.....	42P
<i>Subunit-specific inhibition of recombinantly expressed glycine receptors by ginkgolides and bilobalide</i> R.L. Hawthorne and J.W. Lynch	43P
<i>Crosslinking of $\alpha 1\beta 1$ GABA_A receptor subunits via cysteines introduced into the transmembrane domain</i> T.I. Webb, Z. Yang and J.W. Lynch.....	44P
<i>Etomidate alters the single-channel properties of GABA_A receptors in newborn rat hippocampal neurons</i> V.A.L. Seymour, P.W. Gage and M.L. Tierney	45P
<i>GABA_A $\alpha\beta$ receptors open spontaneously when the conserved M2 leucine 9' residue is mutated to a threonine</i> T.L. Luu, M.L. Tierney and P.W. Gage	46P
<i>GABARAP influences the conductance of recombinant GABA_A channels</i> A.B. Everitt, M.L. Tierney and P.W. Gage.....	47P
<i>C-Terminal peptide of M protein from dengue virus (DVM-C) forms ion channels</i> A. Premkumar, C.R. Horan and P.W. Gage.....	48P
<i>The role of the M1-P1 loop in acid sensitive two-pore domain potassium (TASK) channel regulation</i> C.E. Clarke, A. Mathie and J.I. Vandenberg	49P
<i>Structural studies of chloride intracellular ion channel proteins</i> A.V. Mynott, D.R. Littler, P.M.G. Curmi, L.J. Brown and S.N. Breit	50P
Biomolecular structure and dynamics	
<i>Photochemical behaviour and Na⁺,K⁺-ATPase sensitivity of voltage-sensitive styrylpyridinium fluorescent membrane probes</i> S. Amoroso, V.V. Agon, T. Starke-Peterkovic, M.D. McLeod, H.-J. Apell, P. Sebban and R.J. Clarke	51P
<i>Applications of styrylpyridinium dyes in elucidating ion-transport mechanisms in plant cells</i> S. Amoroso, R.J. Clarke, A. Larkum and R. Quinnell.....	52P
<i>The crystal structure of Pichia pastoris lysyl oxidase at 1.23Å reveals a lysine-lysine covalent cross-link, dehydrolysinonorleucine</i> A.P. Duff, A.E. Cohen, P.J. Ellis, D.B. Langley, D.M. Dooley, H.C. Freeman and J.M. Guss	53P
<i>Molecular dynamics study of conformational changes in human serum albumin by binding of fatty acids</i> S. Fujiwara and T. Amisaki.....	54P
<i>NMR probes of red cell deformation</i> P.W. Kuchel, B.E. Chapman, D.J. Philp and W.A. Bubb	55P
<i>Current-voltage analysis of response to salt stress by salt-tolerant and salt-sensitive charophyte cells</i> M.J. Beilby and V.A. Shepherd.....	56P
<i>Oxygen evolution in chimeric spinach photosystem II with cyanobacteria manganese stabilising protein</i> A. Williamson, W. Hillier, R. Razeghifard and T. Wydrzynski	57P

<i>The role of an oil droplet lens in vision enhancement</i> L. Fischer, M. Vorobyev, A. Zvyagin and T. Plakhotnik	58P
<i>Circular dichroic spectra of the N-terminal region of cardiac myosin binding protein – C</i> C.E. Oakley, L.J. Brown and B.D. Hambly	59P
<i>Changes in mechanical properties of live cell wall during turgor regulation monitored by atomic force microscopy</i> E.M. Mahomudally, M.J. Beilby, V.A. Shepherd and A.R. Moon.....	60P
Cardiac muscle and Skeletal muscle	
<i>Effect of temperature on stretch-induced cardiac action potential shortening in the rat heart: involvement of TREK-1</i> D.R. Kelly, L. Mackenzie and D.A. Saint.....	61P
<i>Immunohistochemical identification of stretch-sensitive two-pore-potassium (TREK) channels in the human heart</i> S.Y. Yuan, H.P. Zhu and D.A. Saint.....	62P
<i>More than one type of stretch activated channel contributes to the action potential duration in guinea pig</i> L. Mackenzie, D.R. Kelly and D.A. Saint.....	63P
<i>Fatty acid composition of red blood cell membranes as a marker of human heart membrane phospholipid fatty acids</i> M.L. Theiss, S. Pepe and P.L. McLennan.....	64P
<i>Confocal Ca²⁺ imaging of mouse sinoatrial node</i> Y.K. Ju, D.G. Allen and M.B. Cannell.....	65P
<i>Reactive oxygen species generated from the mitochondria and not NAD(P)H-oxidase regulate L-type Ca²⁺ channel function during acute hypoxia in ventricular myocytes</i> L.C. Hool, H.M. Viola, C.A. Di Maria and P.G. Arthur	66P
<i>Eccentric damage is accentuated in aged dystrophin-deficient EDL muscles from dystrophic mice (MDX)</i> S. Chan and S.I. Head	67P
<i>Phosphorylation of CSQ affects Ca²⁺ binding and interactions with anchoring protein junctin</i> N.A. Beard, S. Cheung, L. Wei, M. Varsanyi and A.F. Dulhunty.....	68P
<i>A calsequestrin polymer is necessary for the Ca²⁺ binding protein to regulate RyR channels</i> L. Wei, N.A. Beard and A.F. Dulhunty.....	69P
<i>Digoxin and exercise effects on Na⁺,K⁺-pump activity, content, isoform gene and protein expression in human skeletal muscle</i> X. Gong, A.C. Petersen, S. Sostaric, C.A. Goodman, D. Cameron-Smith, R.J. Snow, K.T. Murphy, K. Carey, J. Aw, H. Krum and M.J. McKenna	70P
Nervous system	
<i>Reduced long-term depression is recovered in aging mdx cerebellar Purkinje cells</i> J.L. Anderson, S.I. Head and J.W. Morley.....	71P
<i>Differential action of ω-conotoxins CVID and CVIB on voltage-gated calcium channels in rat sensory neurons</i> L.M. Motin, R.J. Lewis and D.J. Adams	72P

<i>Evidence from collision experiments that onset chopper neurons in the guinea pig cochlear nucleus receive excitatory input from centrifugal collaterals</i> D. Robertson and W.H.A.M. Mulders.....	73P
<i>P2Y receptor activation inhibits the formation and proliferation of primary mouse sub-ventricular-derived neurospheres</i> M.R. Stafford, P.F. Bartlett and D.J. Adams	74P
Membrane transport and its regulation	
<i>Protein kinase A inhibits cell growth induced by overexpression of IK channels</i> C.J. Fowler, K. Ngui, B. Hunne, D. Poole, J.B. Furness and C.B. Neylon.....	75P
<i>Post-transcriptional regulation of CFTR protein expression by 5' untranslated region encoded regulatory elements</i> S-J. Conroy, W.L Davies and A.E.O. Trezise	76P
<i>Nedd4-2, CLC-5 and albumin endocytosis in the proximal tubule: a role for SGK-1?</i> D.H. Hryciw, J. Ekberg, A. Lee, I.L. Lensink, S. Kumar, W.B. Guggino, D.I. Cook, C.A. Pollock and P. Poronnik	77P
<i>Na⁺ H⁺ exchanger regulatory factor 2 (NHERF-2) is a scaffold for the plasma membrane Ca²⁺ ATPase (PMCA)</i> W.A. Kruger, G.R. Monteith, L. Tongpao and P. Poronnik	78P
<i>NHERF1 - a novel scaffold protein for the astroglial glutamate transporter GLAST</i> A. Rayfield, A. Lee, D. Pow, D.H. Hryciw, T.A. Ma, S. Bröer, C.H. Yun and P. Poronnik	79P
<i>Molecular cloning and characterisation of the mouse 'system IMINO' transporter</i> S. Kowalczyk, A. Bröer, M. Munzinger, N. Tietze, K. Klingel and S. Bröer	80P
<i>Increased acetaminophen hepatotoxicity in the NaS1 sulphate transporter null mouse</i> S. Lee, P.A. Dawson and D. Markovich.....	81P
Nutrition	
<i>The influence of dietary fish oil and exercise upon oxidative status biomarkers in a rat model</i> R. Henry, A.J. Owen and P.L. McLennan	82P
Education	
<i>Simulation of visual processing in retinal ganglion cells</i> M. Watson, G. Holmes, T. Byrne and S. Cornford	83P
AuPS Invited Lecture - Angela Dulhunty	
<i>Excitation-contraction coupling from 1969 to 2005</i> A.F. Dulhunty	84P
Symposium 3: Physiology Teaching in the 21st Century: Trends, Challenges and Innovations	
<i>Challenges facing physiology educators in the 21st century</i> Ann Sefton	85P

<i>Research led teaching and learning in physiology</i> R.E. Kemm.....	86P
<i>Problem-based learning (PBL): A novel and effective approach for teaching research skills by addressing contemporary research problems in physiology</i> J. Schwartz and P. Buckley	87P
<i>Improving learning outcomes for students in Clinical Physiology</i> C.R Dallemagne	88P
Symposium 4: Membrane Protein Structure and Interactions	
<i>Seeing spots: miscibility transitions in lipid/cholesterol membranes</i> S.L. Keller and S.L. Veatch.....	89P
<i>Role of the plasma membrane in amyloid formation and toxicity</i> M.I. Aguilar, X. Hou, S. Subasinghe, A. Mechler, K. Hall, L. Martin and D.H. Small	90P
<i>Are chloride intracellular ion channel proteins (CLICs) really channels? Exploring their membrane structure</i> L.J. Brown, D.R. Littler, A.V. Mynott, S.J. Harrop, S.N. Breit, M. Mazzanti and P.M.G. Curmi.....	91P
<i>Structure and dynamics of cellular components using fluorescence and X-ray diffraction techniques</i> Leann Tilley and Nick Klonis	92P
Symposium 5: Regulation of Membrane Transport	
<i>The canonical transient receptor potential channel 1 is an essential structural component of the mechanosensitive calcium permeable channel in vertebrate cells</i> O.P. Hamill, R. Maroto, A. Kurosky, T.G. Wood, A. Raso and B. Martinac	93P
<i>Sensing pressure with K_{2P} channels</i> Eric Honore	94P
<i>Varieties of mechanotransduction: the cytoskeletal stress fibre as a force transmitter and a mechanosensor</i> Masahiro Sokabe, Kimihide Hayakawa and Hitoshi Tatsumi.....	95P
<i>Role of tryptophan residues in ion channel function</i> Amitabha Chattopadhyay	96P
Symposium 6: Membrane Associated Proteins that Regulate Muscle Contraction	
<i>From DHPR to RyR and back again: What lies along the way?</i> Kurt Beam	97P
<i>Regulation of ryanodine receptors from skeletal and cardiac muscle by components of the cytoplasm and lumen</i> D.R. Laver	98P
<i>Structural and functional characterisation of the interaction of the dihydropyridine receptor II-III loop with the ryanodine receptor</i> M.G. Casarotto, Y. Cui, Y. Karunasekara, P.J. Harvey, N. Norris, P.G. Board and A.F. Dulhunty	99P
<i>CLIC-1 chloride channel - matching its properties to a role in skeletal muscle</i> G.Y. Rychkov	100P

Plenary Lecture D. Peter Tieleman

Pores, channels and transporters: computational studies of membrane transport

Peter Tieleman..... 101P

Free Communications 5: Cardiac muscle

Does lignocaine increase the chance of survival from massive heart attack?

S.M. Weiss and P.W. Gage 102P

Acute application of n-3 polyunsaturated fatty acids modify calcium sparks in permeabilised rat cardiac myocytes

B. Honen, R. Dalton, D.F. van Helden and D.R. Laver 103P

Abnormal calcium transients and calcium handling protein expression in cardiomyocytes from mdx (dystrophic) mice

I.A. Williams and D.G. Allen..... 104P

Mechanisms underlying the stretch-dependent slow inotropic response in isolated mouse myocardium

M.L. Ward and D.G. Allen..... 105P

Salutary effects of pyruvate are more evident in female than male glut4-deficient mouse hearts

C.E. Huggins, J.L. Favalaro, J. Proietto, S. Pepe and L.M.D. Delbridge 106P

The angiotensin type 2 receptor prevents cell death in neonatal cardiomyocytes of the hypertrophic heart rat

E.R. Porrello, A. D'Amore, C.L. Curl, S.B. Harrap, W.G. Thomas and L.M.D. Delbridge 107P

Free Communications 6: Cellular signaling

Analyses of the actin cytoskeleton using fluorescence resonance energy transfer (FRET)

C.G. dos Remedios, D. Chhabra, I. Dedova, D. Safer and E DeLaCruz 108P

Profilin binding to sub-micellar concentration of polyphosphoinositides PI(4,5)P₂ and PI(3,4,5)P₃

P.D.J. Moens..... 109P

Phospholipase Cγ is essential for activation of store-operated Ca²⁺ channels in liver cells

T. Litjens, T. Nguyen, E. Aromataris, M. Roberts, G. Barritt and G.Y. Rychkov 110P

Distinct characteristics of exocytosis in mouse pancreatic acinar cells

Peter Thorn, Olga Larina and Ian Parker 111P

Synchronization of Ca²⁺ oscillations through interaction of intracellular Ca²⁺ stores and L-type Ca²⁺ channels

M.S. Imtiaz, J. Zhao, K. Hosaka and D.F. van Helden 112P

Characterization of the electrical activity underlying spontaneous contractions in the mouse ureteropelvic junction

R.J. Lang, B. Zoltkowsky, J. Hammer, W. Meeker, I. Wendt and H.C. Parkington 113P

Free Communications 7: Muscle physiology

Active metabolism of mouse papillary muscle

C. Widén and C.J. Barclay 114P

Functional and electrophoretic identification of two Troponin C isoforms in toad skeletal muscle fibres

B. O'Connell, R. Blazev and G.M.M. Stephenson 115P

<i>X-ray diffraction analysis of the effects of myosin chain-2 phosphorylation on the structure of fast skeletal muscle fibres</i> J.F.Y. Hoh, M. Yamaguchi, M. Kimura, S. Takemori and N. Yagi	116P
<i>Calpain-1 and calpain-3 are not autolysed with exhaustive exercise in humans</i> R.M. Murphy, R.J. Snow and G.D. Lamb	117P
<i>Increased fatigue resistance in EDL muscle of the obese mouse is associated with an increase in the proportion of hybrid IIB+IID fibres</i> R. Blazev, J.G. Kemp, D.G. Stephenson and G.M.M. Stephenson	118P
<i>Insulin-like growth factor-I gene transfer by electroporation enhances skeletal muscle regeneration and function after injury</i> J.D. Schertzer and G.S. Lynch	119P
Free Communications 8: Biophysics	
<i>Bicarbonate is not a physiological substrate of Photosystem II</i> I.L. McConnell, W. Hillier and T. Wydrzynski	120P
<i>Spectra of reef fish – a physics approach to colourful patterns</i> M. Vorobyev and J. Marshall	121P
<i>Resonances of the human vocal tract and some of their uses</i> John Smith and Joe Wolfe	122P
Free communications 9: Education	
<i>Enhancing the first-year experiences of undergraduate students enrolled in large classes</i> R.W. Moni, K.B. Moni, L. Lluca and P. Poronnik	123P
<i>Creating an effective learning community in a large-class service teaching physiology course</i> H.G.G. Ernst and K.L. Colthorpe	124P
<i>The place of physiology in an integrated medical curriculum</i> T.O. Neild	125P
<i>Using a student-centred approach to enhance understanding of the physiology of metabolism and energy balance</i> K.L. Colthorpe and H.G.G. Ernst	126P
<i>Student perceptions and use of pre-specified criteria in constructing complex concept maps in physiology</i> R.W. Moni, E. Beswick, A. Forrest and K.B. Moni	127P
<i>The opinion editorial – a novel assessment task in final year physiology</i> D.H. Hryciw, P. Poronnik and R.W. Moni	128P
Free communications 10: Skeletal Muscle 1	
<i>Low dose formoterol treatment reverses sarcopenia and improves muscle function in fast- but not slow-twitch skeletal muscles of aged rats</i> G.S. Lynch and J.G. Ryall	129P
<i>β-adrenergic signalling in skeletal muscle regeneration after myotoxic injury</i> F. Beitzel, M.N. Sillence and G.S. Lynch	130P

<i>Streptomycin reduces stretch-induced membrane permeability in isolated muscles from mdx (dystrophic) mice</i> N.P. Whitehead, M. Streamer and D.G. Allen	131P
<i>Muscle weakness in a mouse model of nemaline myopathy can be reversed with exercise and reveals a novel myofibre repair mechanism</i> A.J. Kee, J.E. Joya, V. Nair-Shalliker, M.-A. Nguyen, M. Ghodduzi and E.C. Hardeman.....	132P
<i>Contraction-mediated damage in mdx dystrophic mouse tibialis anterior muscles is not affected by the membrane sealant poloxamer</i> D.R. Plant, J.G. Ryall and G.S. Lynch.....	133P
<i>Aberrant splicing of ryanodine receptor reduces Ca²⁺ release via an inter-domain interaction in myotonic dystrophy type 1</i> T. Kimura, M. Nakamori, J.D. Lueck, P. Pouliquin, R.T. Dirksen, M.P. Takahashi, S. Sakoda and A.F. Dulhunty.....	134P
Free communications 11: Skeletal muscle 2	
<i>The effect of altering the rest period during interval training on adaptations to muscle metabolism, ion regulation and exercise performance</i> J. Edge, D. Bishop and C.A. Goodman.....	135P
<i>Long lasting muscle fatigue: partial disruption of EC-coupling by the elevated cytosolic calcium during contractions</i> E. Verburg, T.L. Dutka and G.D. Lamb	136P
<i>The role of reactive oxygen species on stretch-induced muscle damage in dystrophic mice</i> D.G. Allen and E.Y. Yeung	137P
<i>Effects of raising the temperature from 25°C to 37°C on twitch responses in fast-twitch mechanically skinned muscle fibres of the rat</i> C. van der Poel, J. Edwards and D.G. Stephenson.....	138P
<i>Exposure of mammalian skeletal muscle to sub-physiological temperatures reduces its ability to function at physiological temperatures</i> J. Edwards, C. van der Poel and D.G. Stephenson.....	139P
Free communications 12: Systems physiology	
<i>Novel nifedipine-insensitive high voltage activated calcium channels play a role in vascular tone of cerebral arteries</i> M.F. Navarro-Gonzalez and C.E. Hill.....	140P
<i>Accentuation during diabetes of differential connexin expression between the preglomerular and postglomerular renal vasculature</i> J.H. Zhang and C.E. Hill.....	141P
<i>High-amplitude oscillations in human skin blood flow are distinct from known cardiac or respiratory influences</i> C.D. Haley, C.J. Gordon, N.A.S. Taylor and A.B. Jenkins.....	142P
<i>Activation of at least three classes of ion channels by β-adrenoceptor activation in pregnant uterine smooth muscle</i> H.C. Parkington, M.A. Tonta, S. Simon, S.A. Cohen, A. Satragno, R.J. Lang and H.A. Coleman	143P

<i>Expression of a constitutively active K⁺ channel prevents cell division in the mouse preimplantation embryo</i> M.L. Day, C.G. Bailey, J.E. Rasko, M.H. Johnson and D.I. Cook	144P
--	------

Symposium 7: Function and Regulation of Ion Transport Membrane Proteins

<i>The transition between gating states in inward rectifier K⁺ channels</i> J. Gulbis	145P
---	------

<i>Channelrhodopsin 1,2, a new class of ion channels: functional description and cellular applications</i> Georg Nagel, Peter Hegemann, Suneel Kateriya and Ernst Bamberg	146P
--	------

<i>Membrane lipid composition and its effect on Na⁺,K⁺-ATPase molecular activity: insights from mammals, birds and ectotherms</i> N. Turner, P.L. Else, B.J. Wu and A.J. Hulbert	147P
---	------

<i>Testing the membrane pacemaker model of metabolism</i> P.L. Else, N. Turner, T.W. Mitchell, B.J. Wu and A.J. Hulbert	148P
--	------

<i>Regulation of the Na,K-ATPase</i> Helge Rasmussen	149P
---	------

Symposium 8: Epithelial Transport of Ions and Metabolites

<i>A systems biology approach to understanding the role of peptide transporters in biology</i> H. Daniel, B. Meissner, B. Spanier, D. Weitz and I. Frey	150P
--	------

<i>Na⁺-H⁺ exchange regulatory factors NHERF-1 and NHERF-2: roles in albumin endocytosis in the proximal tubule</i> P. Poronnik, C. Ferguson, R. Parton, C.H. Yun and D.H. Hryciw	151P
---	------

<i>Sulphate ions in mammalian physiology: lessons from sulphate transporter knock-out mice</i> P.A. Dawson, B. Gardiner, S. Lee, M.C. Ku, S.M. Grimmond and D. Markovich	152P
---	------

<i>Regulation of epithelial Na⁺ channels</i> A. Dinudom, I.-H. Lee and D.I. Cook	153P
--	------

<i>Disorders of neutral amino acid resorption in epithelial cells</i> S. Bröer	154P
---	------